

# ABHD14A Antibody



PACO07529

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ABHD14A

**Uniprot**

Q9BUJ0

**Synonyms:**

abhydrolase domain containing 14A; ABHD14A; DKFZp564O243; DORZ1; FLJ60317; FLJ60467 ;

**Immunogen:**

Human ABHD14A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A