

# ABHD2 Antibody



PACO07531

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ABHD2

**Uniprot**

P08910

**Synonyms:**

abhydrolase domain containing 2; ABHD2; HS1-2; LABH2; MGC111112; MGC26249; PHP51-2 ;

**Immunogen:**

Human ABHD2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A