

ABTB1 Antibody



PACO07538

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ABTB1

Uniprot

Q969K4

Synonyms:

ankyrin repeat and BTB (POZ) domain containing 1; ABTB1; BPOZ; Btb3; EF1ABP; MGC20585; PP2259 ;

Immunogen:

Human ABTB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A