

ACADSB Antibody



PACO07546

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ACADSB

Uniprot

P45954

Synonyms:

acyl-Coenzyme A dehydrogenase, short/branched chain; ACADSB;2-MEBCAD; ACAD7; SBCAD ;

Immunogen:

Human ACADSB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A