

ACADVL Antibody



PACO07547

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ACADVL

Uniprot

P49748

Synonyms:

acyl-Coenzyme A dehydrogenase, very long chain; ACADVL; ACAD6; LCACD; VLCAD ;

Immunogen:

Human ACADVL

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A