

ACAP2 Antibody



PACO07549

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ACAP2

Uniprot

Q15057

Synonyms:

ArfGAP with coiled-coil, ankyrin repeat and PH domains 2; ACAP2; CENTB2; CNT-B2; KIAA0041 ;

Immunogen:

Human ACAP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A