

ACAP3 Antibody



PACO07550

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ACAP3

Uniprot

Q96P50

Synonyms:

ArfGAP with coiled-coil, ankyrin repeat and PH domains 3; ACAP3; CENTB5; KIAA1716 ;

Immunogen:

Human ACAP3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A