

ASIC4 Antibody



PACO07554

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ASIC4

Uniprot

Q96FT7

Synonyms:

amiloride-sensitive cation channel 4, pituitary; ACCN4; ASIC4; BNAC4; MGC17248; MGC24860 ;

Immunogen:

Human ACCN4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A