

ACOT11 Antibody



PACO07561

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ACOT11

Uniprot

Q8WXI4

Synonyms:

acyl-CoA thioesterase 11; ACOT11; BFIT; BFIT1; BFIT2; KIAA0707; STARD14; THEA; THEM1 ;

Immunogen:

Human ACOT11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A