

ACOT7 Antibody



PACO07563

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ACOT7

Uniprot

O00154

Synonyms:

acyl-CoA thioesterase 7; ACOT7; ACH1; ACT; BACH; CTE-II; LACH; LACH1; MGC1126; hBACH ;

Immunogen:

Human ACOT7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A