

ACSS1 Antibody



PACO07575

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ACSS1

Uniprot

Q9NUB1

Synonyms:

acyl-CoA synthetase short-chain family member 1; ACSS1; ACAS2L; AceCS2L; FLJ45659; MGC33843 ;

Immunogen:

Human ACSS1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A