

ACTN2 Antibody



PACO07585

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ACTN2

Uniprot

P35609

Synonyms:

actinin, alpha 2; ACTN2; CMD1AA ;

Immunogen:

Human ACTN2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A