

ACTR5 Antibody



PACO07593

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ACTR5

Uniprot

Q9H9F9

Synonyms:

ARP5 actin-related protein 5 homolog (yeast); ACTR5; Arp5; FLJ12785; INO80M ;

Immunogen:

Human ACTR5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A