

ACYP1 Antibody



PACO07598

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ACYP1

Uniprot

P07311

Synonyms:

acylphosphatase 1, erythrocyte (common) type; ACYP1; ACYPE ;

Immunogen:

Human ACYP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A