

ACYP2 Antibody



PACO07599

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ACYP2

Uniprot

P14621

Synonyms:

acylphosphatase 2, muscle type; ACYP2; ACYM; ACYP ;

Immunogen:

Human ACYP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A