

ADAD2 Antibody



PACO07600

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ADAD2

Uniprot

Q8NCV1

Synonyms:

adenosine deaminase domain containing 2; ADAD2; FLJ00337; TENRL ;

Immunogen:

Human ADAD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A