

ADAM12 Antibody



PACO07601

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ADAM12

Uniprot

O43184

Synonyms:

ADAM metalloproteinase domain 12; ADAM12; MCMP; MCMPMItna; MLTN; MLTNA ;

Immunogen:

Human ADAM12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A