

ADAT2 Antibody



PACO07606

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ADAT2

Uniprot

Q7Z6V5

Synonyms:

adenosine deaminase, tRNA-specific 2, TAD2 homolog (S. cerevisiae); ADAT2; DEADC1; DKFZp686L1118; FLJ44213; TAD2; dJ20N2; dJ20N2.1 ;

Immunogen:

Human ADAT2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A