

ADCY6 Antibody



PACO07610

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ADCY6

Uniprot

O43306

Synonyms:

adenylate cyclase 6; ADCY6; AC6; DKFZp779F075; KIAA0422 ;

Immunogen:

Human ADCY6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A