

ADD1 Antibody



PACO07612

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ADD1

Uniprot

P35611

Synonyms:

adducin 1 (alpha); ADD1; ADDA; MGC3339; MGC44427 ;

Immunogen:

Human ADD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A