

ADH6 Antibody



PACO07617

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ADH6

Uniprot

P28332

Synonyms:

alcohol dehydrogenase 6 (class V); ADH6; ADH-5 ;

Immunogen:

Human ADH6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A