

# ADHFE1 Antibody



PACO07618

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ADHFE1

**Uniprot**

Q8IWW8

**Synonyms:**

alcohol dehydrogenase, iron containing, 1; ADHFE1; ADH8; FLJ32430; HMFT2263; HOT; MGC48605 ;

**Immunogen:**

Human ADHFE1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A