

ADRM1 Antibody



PACO07629

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ADRM1

Uniprot

Q16186

Synonyms:

adhesion regulating molecule 1; ADRM1; ARM1; GP110; MGC29536; Rpn13 ;

Immunogen:

Human ADRM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A