

AFAP1 Antibody



PACO07633

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AFAP1

Uniprot

Q8N556

Synonyms:

actin filament associated protein 1; AFAP1; AFAP; AFAP-110; FLJ56849 ;

Immunogen:

Human AFAP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A