

AFG3L2 Antibody



PACO07636

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

AFG3L2

Uniprot

Q9Y4W6

Synonyms:

AFG3 ATPase family gene 3-like 2 (yeast); AFG3L2; FLJ25993 ;

Immunogen:

Human AFG3L2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A