

AFP Antibody



PACO07637

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AFP

Uniprot

P02771

Synonyms:

alpha-fetoprotein; AFP; FETA; HPAFP ;

Immunogen:

Human AFP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A