

AGBL2 Antibody



PACO07639

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

AGBL2

Uniprot

Q5U5Z8

Synonyms:

ATP/GTP binding protein-like 2; AGBL2; FLJ23598 ;

Immunogen:

Human AGBL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A