

AGPAT6 Antibody



PACO07644

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GPAT4

Uniprot

Q86UL3

Synonyms:

1-acylglycerol-3-phosphate O-acyltransferase 6 (lysophosphatidic acid, acyltransferase, zeta); AGPAT6; DKFZp586M1819; LPAAT-zeta; LPAATZ; TSARG7 ;

Immunogen:

Human AGPAT6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A