

PACO07648

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

AGT

**Uniprot**

P01019

**Synonyms:**

angiotensinogen (serpin peptidase inhibitor, clade A, member 8); AGT; ANHU; FLJ92595; FLJ97926; SERPINA8 ;

**Immunogen:**

Human AGT.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A