

AGTPBP1 Antibody



PACO07649

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AGTPBP1

Uniprot

Q9UPW5

Synonyms:

ATP/GTP binding protein 1; AGTPBP1; DKFZp686M20191; KIAA1035; NNA1 ;

Immunogen:

Human AGTPBP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A