

# AK3 Antibody



PACO07662

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

AK3

**Uniprot**

Q9UIJ7

**Synonyms:**

adenylate kinase 3; AK3; AK3L1; AK6; AKL3L; AKL3L1; FIX ;

**Immunogen:**

Human AK3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A