

AKR1A1 Antibody



PACO07670

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

AKR1A1

Uniprot

P14550

Synonyms:

aldo-keto reductase family 1, member A1 (aldehyde reductase); AKR1A1; ALDR1; ALR; ARM; DD3; MGC12529; MGC1380 ;

Immunogen:

Human AKR1A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A