

AKR1C3 Antibody



PACO07672

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AKR1C3

Uniprot

P42330

Synonyms:

aldo-keto reductase family 1, member C3 (3-alpha hydroxysteroid dehydrogenase, type II); AKR1C3; DD3; DDX; HA1753; HAKRB; HAKRe; HSD17B5; KIAA0119; hluPGFS ;

Immunogen:

Human AKR1C3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A