

ALDH18A1 Antibody



PACO07679

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ALDH18A1

Uniprot

P54886

Synonyms:

aldehyde dehydrogenase 18 family, member A1; ALDH18A1; GSAS; MGC117316; P5CS; PYCS ;

Immunogen:

Human ALDH18A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A