

# ALDH1B1 Antibody



PACO07682

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ALDH1B1

**Uniprot**

P30837

**Synonyms:**

aldehyde dehydrogenase 1 family, member B1; ALDH1B1; ALDH5; ALDHX; MGC2230 ;

**Immunogen:**

Human ALDH1B1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A