

# ALDH8A1 Antibody



PACO07691

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ALDH8A1

**Uniprot**

Q9H2A2

**Synonyms:**

aldehyde dehydrogenase 8 family, member A1; ALDH8A1; ALDH12; DJ352A20.2; DKFZp779D2315; MGC138650 ;

**Immunogen:**

Human ALDH8A1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A