

# ALG1 Antibody



PACO07695

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ALG1

**Uniprot**

Q9BT22

**Synonyms:**

asparagine-linked glycosylation 1, beta-1,4-mannosyltransferase homolog (S. cerevisiae); ALG1; HMAT1; HMT-1; HMT1 ;

**Immunogen:**

Human ALG1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A