

ALG13 Antibody



PACO07696

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

ALG13

Uniprot

Q9NP73

Synonyms:

asparagine-linked glycosylation 13 homolog (S. cerevisiae); ALG13; CXorf45; FLJ23018; GLT28D1; MDS031; MGC12423; YGL047W ;

Immunogen:

Human ALG13.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A