

ALG6 Antibody



PACO07698

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ALG6

Uniprot

Q9Y672

Synonyms:

asparagine-linked glycosylation 6, alpha-1,3-glucosyltransferase homolog (S. cerevisiae);

Immunogen:

Human ALG6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A