

# AMICA1 Antibody



PACO07714

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

JAML

**Uniprot**

Q86YT9

**Synonyms:**

adhesion molecule, interacts with CXADR antigen 1; AMICA1; AMICA; CREA7-1; CREA7-4; FLJ37080; Gm638; JAML; MGC118814; MGC118815 ;

**Immunogen:**

Human AMICA1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A