

AMIGO2 Antibody



PACO07715

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AMIGO2

Uniprot

Q86SJ2

Synonyms:

adhesion molecule with Ig-like domain 2; AMIGO2; ALI1; DEGA ;

Immunogen:

Human AMIGO2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A