

AMN1 Antibody



PACO07717

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AMN1

Uniprot

Q8IY45

Synonyms:

antagonist of mitotic exit network 1 homolog (S. cerevisiae);

Immunogen:

Human AMN1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A