

# ANAPC4 Antibody



PACO07728

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ANAPC4

**Uniprot**

Q9UJX5

**Synonyms:**

anaphase promoting complex subunit 4; ANAPC4; APC4 ;

**Immunogen:**

Human ANAPC4.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A