

ANGPTL3 Antibody



PACO07733

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ANGPTL3

Uniprot

Q9Y5C1

Synonyms:

angiopoietin-like 3; ANGPTL3; ANGPT5 ;

Immunogen:

Human ANGPTL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A