

ANGPTL7 Antibody



PACO07736

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ANGPTL7

Uniprot

O43827

Synonyms:

angiopoietin-like 7; ANGPTL7; AngX; CDT6; RP4-647M16.2; dJ647M16.1 ;

Immunogen:

Human ANGPTL7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A