

ANKMY2 Antibody



PACO07737

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ANKMY2

Uniprot

Q8IV38

Synonyms:

ankyrin repeat and MYND domain containing 2; ANKMY2; DKFZp564O043; ZMYND20 ;

Immunogen:

Human ANKMY2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A