

ANO3 Antibody



PACO07742

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ANO3

Uniprot

Q9BYT9

Synonyms:

anoctamin 3; ANO3; C11orf25; GENX-3947; TMEM16C ;

Immunogen:

Human ANO3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A