

ANTXR1 Antibody



PACO07747

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ANTXR1

Uniprot

Q9H6X2

Synonyms:

anthrax toxin receptor 1; ANTXR1; ATR; FLJ10601; FLJ11298; FLJ21776; TEM8 ;

Immunogen:

Human ANTXR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A