

ANXA8L2 Antibody



PACO07758

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ANXA8L1

Uniprot

Q5VT79

Synonyms:

annexin A8-like 2; ANXA8L2; ANXA8; FLJ32754; FLJ54151; bA145E20.2 ;

Immunogen:

Human ANXA8L2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A