

ANXA9 Antibody



PACO07759

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ANXA9

Uniprot

O76027

Synonyms:

annexin A9; ANXA9; ANX31 ;

Immunogen:

Human ANXA9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A